

DispensMate plus Bottle-top Dispenser Chemical Compatibility at 20°C

The devices of SCIOLOGEX-DispensMate plus which contact with dispensed liquid consist of BSG, PTFE, FEP, and closure cap of outlet is PP; non contact liquids parts consist of PC and other materials. Please note that these tables are just a directional guide not the manufacturer's commitment. Please read the user manual carefully before use and to do related experiments necessarily, which can be used to determine whether should be used. Good laboratory practice would be to rinse out the liquid handing unit at the end of each day with distilled water to prevent corrosive liquids being left in contact with the parts for too long.

We referred to the general technical data and public information from related companies. The table below is not our proprietary data, for user's reference only.

CHEMICAL	BSG	PTFE	FEP	PC	PP
Acids					
Acetic, Glacial	R			NR	R
Acetic, 25%	R	R	R	R	R
Hydrochloric, Concentrated	R				
Hydrochloric, 20%	R	R	R	SR	R
Sulphuric, concentrated	R				
Sulphuric, 25%	R	R	R	R	R
Nitric, Concentrated	R				
Nitric, 30%	R	R	R	R	SR
Phosphoric, 25%	R	R	R		4
Formic, 25%	R	R	R		
Trichloroacetic, 10%	R	R	R	SR	SR
Formic Acid, 85%	R	R	R	R	R
Arsenic Acid	R				
Boric Acid, 10%	R	R	R	R	R
Chromic Acid, 10%	R	R	R	NR	R
Hydrofluoric Acid, 35%	NR	Exceptions	R	NR	R
Phosphoric Acid 85%	R	R	R	R	R
Nitric Acid, 50%	R	R	R		
Sulphuric Acid, 95%	R	R	R	NR	NR
Alkalies					
AmmoniumHydroxide,25%	R	R	R	NR	R
Potassium Hydroxide	R	R	R	NR	R
Sodium Hydroxide	R	R	R	NR	R
Alcohols					
Methanol, 98%	R	R		R	R
Ethanol, 98%	R			R	R
Ethanol, 70%	R			R	R
Isopropanol,n-Propanol	R			R	R
Amyl Alcohol, Butanol	R				
Benzyl Alcohol	R	R	R	SR	SR
Ethylene Glycol	R	R	R	R	R
Propylene Glycol	R	R	R	R	R
Glycerol	R	R	R	R	R
Hydrocarbons					
Hexane, Xylene	R	R	R	NR	R
Toluene, Benzene	R	R	R	NR	SR
Kerosene, Gasoline	R				
Tetralin, Decalin	R				
Halogenated Hydrocarbons					
Methyl Chloride	R			NR	SR
Chloroform	R	R	R	NR	NR
Trichloroethylene	R	R	R	NR	NR
Monochlorobenzene, Freon	R				
Carbon Tetrachloride	R	R	R	NR	NR

Ketones						
Acetone	R	R	R	NR	R	
Methyl Ethyl Ketone	R	R				
Isopropylacetone	R					
Methyl Isobutyl Ketone	R					
Esters						
Ethyl Acetate	R	R		NR	R	
Methyl Acetate	R					
Amyl & Propyl Acetate	R					
Butyl Acetate	R	R	R	NR	NR	
Propylene Glycol Acetate	R					
2-Ethoxyethyl Acetate	R					
Methyl Cellosolve Acetate	R					
Benzyl Benzoate	R					
Isopropyl Myristate	R					
Tricresol Phosphate	R					
Oxides – Ethers						
Ethyl Ether	R					
1,4 Dioxane & Tetrahydrofuran	R	R	R	NR	SR	
Dimethylsulphoxide (DMSO)	R	R	R	NR	R	
Isopropyl Ether	R			NR	NR	
Solvents with Nitrogen						
Dimethyl Formamide	R	R	R			
Diethylacetamide	R	R				
Triethanolamine	R					
Aniline	R	R	R	SR	R	
Pyridine	R	R	R	NR	SR	
Miscellaneous						
Phenol, Aqueous, 10%	R					
Formaldehyde Solution, 30%	R	R	R	R	R	
Hydrogen Peroxide, 30%	R	R	R	R	R	
Silicone Oil & Mineral Oil	R					
Pyridine	R	R	R	NR	SR	
Acetaldehyde	R	R	R	SR	R	
Ammonia, 25% ac. Sol.	R	R		NR	R	
Ammonium	R					
Calcium Chloride aq. Sol.	R	R	R	R	R	
Chlorine	R	R	R			
Chlorobenzene	R			NR	NR	
Fluorinated Hydrocarbons	R					
Hexane	R	R	R	R	R	
Iodine (tincture of)	R	R				
Potassium Chloride aq. Sol.	R			R	R	
Potassium Permanganate aq. Sol.	R			R	R	
Magnesium Chloride aq. Sol.	R					
Methylene Chloride	R	R	R	NR	SR	
Sodium Carbonate	R					
Sodium Dichromate	R	R	R	R	R	
Phenol, 100%	R	R	R	NR	R	
Mercury	R	R	R	R	R	
Silver Nitrate	R	R	R	R	R	
Toluene	R	R	R	NR	SR	
Hydrogen Peroxide, 30%	R	R	R	NR	R	
Xylene	R	R	R	NR	NR	
Zinc Chloride, 10%	R	R	R	R	R	
Zinc Sulphate, 10%	R	R	R	R	R	

KEY:
R = RESISTANT **NR** = NON-RESISTANT **SR** = SLIGHTLY RESISTANT
EXCEPTIONS = RESISTANT WITH EXCEPTIONS