
# Kühlwasserregelventil / Cooling water control valve #10312

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## ENGLISH

The cooling water control valve controls the cooling water consumption according to the demand and, consequently, reduces it. The Pilot ONE controller controls the cooling water control valve so that only the cooling water volume flows through the **>cooling coil< [29]** that is currently required for the temperature control task.

See also Section »**Recommended temperature control and cooling water hoses**« in the operation manual of the employed temperature control unit.

Permissible temperature range of the cooling water	-20 °C bis +80 °C
Permissible maximum temperature of the bath	≤ 95 °C
Cooling water medium See also Section » <b>Temperature control units with water cooling</b> « in the operation manual of the employed temperature control unit.	a) Water (municipal water or filtered fluvial water) b) Water-MEG mixture
Cooling water hoses	Designed for mains water pressure
Maximum water pressure	6 bars
Protection rating EN 60529	IP20

## INFORMATION

In the event of a mains failure, the cooling water control valve no longer controls the cooling water consumption. The cooling water control valve remains in its last position (open or closed).

The cooling water consumption is no longer controlled if the cooling water control valve is **open** when power fails. Ensure that the cooling water return line can accommodate the water volume of a completely open cooling water control valve.

If the cooling water control valve is **closed** when power fails, no cooling water is routed through the **>cooling coil< [29]**.

## Installing the cooling water control valve

### PROCEDURE

#### Securing the cooling water control valve

- Ensure that the temperature control unit is switched off and de-energized.
- Install the cooling water control valve on a wall in the vicinity of the temperature control unit. Use the **>mounting opening< [C]** in the enclosure of the cooling water control valve for this purpose.

### PROCEDURE

#### Connecting the cooling water control valve with the >cooling coil< [29]

- Connect the **>cooling water connection< [A]** at the cooling water control valve with the inlet of the **>cooling coil< [29]** you want to use. Observe the specified **>installation direction< [E]** for the connection. Do not install the cooling water control valve at the **>cooling coil< [29]** outlet you use.
- Connect the **>cooling water connection< [B]** at the cooling water control valve with the cooling water supply of the building.
- Check the connections for leaks.

### PROCEDURE

#### Connecting the cooling water control valve with the temperature control unit

- Ensure that the temperature control unit is switched off and de-energized.
- Connect the **>control cable< [D]** on the cooling water control valve with the **>service interface< [50]** connection on the temperature control unit.

## Switching on the temperature control unit

### PROCEDURE

- Switch on the temperature control unit at the **>mains isolator< [36]/>mains switch< [37]**. Provided that the cooling water control valve was connected correctly, the **>touchscreen< [88]** shows the message "Cooling water control valve active" in the **>status cell< [Field 10]**.