



CUSTOM FLASKS

CUSTOM FABRICATED HEAVY WALL FLASKS


Chemglass can custom fabricate round bottom, flat bottom, pear shaped, recovery, european, and erlenmeyer flasks using our exclusive HEAVY WALL flask blanks. Flasks can be supplied with any number of the special features listed below. Please contact our Technical Service Department at 800.843.1794 or via e-mail at technical-service@chemglass.com with the following information, and we will provide you with a quotation. A sketch of your design can also be faxed to 856.696.7366.




CG-1539




CG-1538-M





Location of Side Necks



CG-1968-Q Drain



Stopcock Drain



Flask Type:

Round Bottom Recovery
 Flat Bottom European
 Pear Shaped Erlenmeyer
 Other: _____

Capacity:

25mL	250mL	6,000mL
50mL	500mL	12,000mL
100mL	1,000mL	22,000mL
125mL	2,000mL	50,000mL
150mL	3,000mL	72,000mL
200mL	5,000mL	100,000mL

Neck Options:

Outer Joint O-Ring Joint
 Socket Joint GL Thread
 Internal Thread Other: _____
 Schott O-Ring Flange

For A Complete Listing Of Joint Types And Sizes, See Pages 265-290.

Center Neck:

Size: 24/40, 65/40, etc...
 Type: Outer Joint, Socket Joint, etc...

Side Neck(s):

How Many: 1, 2, 3, etc...
 Size: 24/40, 65/40, etc...
 Type: Outer Joint, Socket Joint, etc...
 Location: 3:00, 6:00, 9:00, Or 12:00
 Straight or Angled: Standard Angle Is 20°

Bottom Outlet:

CG-510 Chem-Drain Valve (6mm or 12mm)
 CG-511 Chem-Cap Valve (12mm or 20mm)
 CG-1968-Q Beaded Pipe with Detachable Valve
 Teflon® Stopcock (2mm, 4mm, etc...)
 Glass Stopcock (2mm, 4mm, etc...)

Flask Options:

Plastic Coating (Temperature Range -20 To 140°F)
 Morton Indents (Not Suitable For Vacuum)
 Half Jacket **
 Full Jacket **

** For Safety Purposes, GL Threads With PP Or Teflon® Detachable Hose Connectors Can Be Used On Inlet/Outlet Instead Of Standard Glass Serrated Hose Connections

Special Instructions:

Quantity: