



CG-1605-A FLOWMETER, NITROGEN, METERING VALVE

Precision machined clear acrylic Nitrogen flowmeter with Brass metering valve is for use with Chemglass process reactors. Front of flowmeter has a blue scale indicating the L/min of air. HDPE hose connection is for use with 5/16" ID tubing. Unit is supplied complete with a 1/2" OD mounting arm for ease in clamping to support frame. Flowmeters with larger L/min ranges are available. Please contact us for details.

Specifications:

Hose Barb OD (mm):	8
Dimensions W x D x H (mm):	25.4 x 28.6 x 102
MAX Operating Pressure (PSIG/BAR):	100/6.89
MAX Operating Temperature (°C/°F):	65/150
Body Material:	Clear Acrylic
Valve Material:	Brass
Seals Material:	Buna-N
Accuracy:	± 5% Full Scale

Part Number	Range (L/min) Air
CG-1605-A-10	0.04 to 0.5
CG-1605-A-12	0.4 to 5.0
CG-1605-A-14	1 to 10.0



CG-1605-A



Custom Glassware

Chemglass has the ability to custom fabricate glassware to your unique specifications. We specialize in custom:

- Chromatography Columns
- Flasks
- Manifolds
- Condensers
- Distillation apparatus
- Funnels
- Reaction Vessels
- and many more...

Contact our Technical Service Department at 01525 382277 or email: info@gpelimited.co.uk.