

CG-100 GROUND JOINT, INNER, FULL LENGTH

Standard taper inner joint, full length.

Part Number	Joint Size	Tubing OD x Wall (mm)
CG-100-01	7/25	6.35 x 1.2
CG-100-02	10/30	8.00 x 1.5
CG-100-03	12/30	10.0 x 1.5
CG-100-04	14/35	12.7 x 1.6
CG-100-05	19/38	17.0 x 1.8
CG-100-06	24/40	22 x 2.1
CG-100-07	29/42	28 x 2.0
CG-100-08	34/45	32 x 2.4
CG-100-09	40/50	38 x 2.0
CG-100-10	45/50	41 x 2.0
CG-100-11	50/50	48 x 2.0
CG-100-12	55/50	51 x 2.0
CG-100-13	60/50	57 x 2.4
CG-100-14	71/60	70 x 2.4
CG-100-15	103/60	100 x 2.4




CG-100

CG-101 GROUND JOINT, OUTER, FULL LENGTH


Standard taper outer joint tooled on tubing approximately the same diameter as the top reinforcing ring, affording no hold-up below joint.

Part Number	Joint Size	Tubing OD x Wall (mm)
CG-101-01	7/25	10 x 1.5
CG-101-02	10/30	14 x 1.2
CG-101-03	12/30	16 x 1.2
CG-101-04	14/35	19 x 1.8
CG-101-05	19/38	22 x 1.8
CG-101-06	24/40	28 x 2.5
CG-101-07	29/42	32 x 2.4
CG-101-08	34/45	38 x 2.0
CG-101-09	40/50	45 x 2.0
CG-101-10	45/50	48 x 2.3
CG-101-11	50/50	54 x 2.4
CG-101-12	55/50	60 x 2.4
CG-101-13	60/50	64 x 2.4
CG-101-14	71/60	75 x 2.4
CG-101-15	103/60	110 x 3.0



CG-101

JOINTS, GROUND



CG-102

CG-102 GROUND JOINT, OUTER, FULL LENGTH

Standard taper outer joint tooled on the same size tubing as the CG-100 inner joints. All CG-102 outer joints are supplied standard with rings, but are available tooled without rings on a special order basis.

Part Number	Joint Size	Tubing OD x Wall (mm)
CG-102-01	7/25	6.35 x 1.2
CG-102-02	10/30	8.00 x 1.5
CG-102-03	12/30	10.0 x 1.5
CG-102-04	14/35	12.7 x 1.6
CG-102-05	19/38	17.0 x 1.8
CG-102-06	24/40	22 x 1.8
CG-102-07	29/42	28 x 2.0
CG-102-08	34/45	32 x 2.4
CG-102-09	40/50	38 x 2.0
CG-102-10	45/50	41 x 2.0
CG-102-11	50/50	48 x 2.0
CG-102-12	55/50	51 x 2.0
CG-102-13	60/50	57 x 2.4
CG-102-14	71/60	70 x 2.4



CG-103

CG-103 GROUND JOINT, OUTER, FULL LENGTH

Standard taper outer joints tooled with reinforcing rings.

Part Number	Joint Size	Tubing OD x Wall (mm)
CG-103-01	7/25	10 x 1.5
CG-103-02	10/30	14 x 1.2
CG-103-03	12/30	16 x 1.2
CG-103-04	14/35	19 x 1.8
CG-103-05	19/38	22 x 1.8
CG-103-06	24/40	28 x 2.5
CG-103-07	29/42	32 x 2.4
CG-103-08	34/45	38 x 2.0
CG-103-09	40/50	45 x 2.0
CG-103-10	45/50	48 x 2.0
CG-103-11	50/50	54 x 2.4
CG-103-12	55/50	60 x 2.4
CG-103-13	60/50	64 x 2.4
CG-103-14	71/60	75 x 2.4

CG-104 GROUND JOINT, INNER, MEDIUM LENGTH

Part Number	Joint Size	Tubing OD x Wall (mm)
CG-104-01	5/5	4.00 x 0.8
CG-104-02	5/12	4.00 x 0.8
CG-104-03	7/10	6.35 x 1.2
CG-104-04	7/15	6.35 x 1.2
CG-104-05	10/18	8.00 x 1.5
CG-104-06	12/18	10.0 x 1.5
CG-104-07	14/10	12.7 x 1.6
CG-104-08	14/20	12.7 x 1.6
CG-104-09	19/22	17.0 x 1.8
CG-104-10	24/25	22.0 x 1.8
CG-104-11	29/26	28 x 2.0
CG-104-12	34/28	32 x 2.4
CG-104-13	40/35	38 x 2.0



CG-104

CG-105 GROUND JOINT, OUTER, MEDIUM LENGTH

Part Number	Joint Size	Tubing OD x Wall (mm)
CG-105-01	5/12	8.0 x 1.5
CG-105-02	7/10	10 x 1.5
CG-105-03	7/15	10 x 1.5
CG-105-04	10/18	14 x 1.2
CG-105-05	12/18	16 x 1.2
CG-105-06	14/20	19 x 1.8
CG-105-07	19/22	22 x 1.8
CG-105-08	24/25	28 x 2.0
CG-105-09	29/26	32 x 2.4
CG-105-10	34/28	38 x 2.0
CG-105-11	40/35	45 x 2.0



CG-105

CG-106 GROUND JOINT, OUTER, MEDIUM LENGTH

Standard taper outer joint tooled on the same diameter tubing as the CG-104 inner joints.

Part Number	Joint Size	Tubing OD x Wall (mm)
CG-106-01	5/12	7.0 x 1.0
CG-106-02	7/15	7.0 x 1.0
CG-106-03	10/18	8.0 x 1.5
CG-106-04	12/18	10.0 x 1.5
CG-106-05	14/20	12.7 x 1.6
CG-106-06	19/22	17 x 1.8
CG-106-07	24/25	22 x 1.8
CG-106-08	29/26	28 x 2.0
CG-106-09	34/28	32 x 2.4
CG-106-10	40/35	38 x 2.0



CG-106