

MELTING POINT TUBES

CG-1840 MELTING POINT TUBE, OPEN BOTH ENDS

Tube is manufactured from borosilicate glass and packaged 100 tubes per vial. The Code -02 has a smaller ID to conform to USP XXI requirements.

Part Number	Tube ID (mm)	Wall (mm)	Length (mm)	Vial Qty
CG-1840-01	1.5-1.8	0.20	100	100
CG-1840-02	0.8-1.1	0.25	100	100

CG-1841 MELTING POINT TUBE, CLOSED ONE END


Similar to CG-1840, but closed at one end. Tube is manufactured from borosilicate glass and packaged 100 tubes per vial. The Code -02 has a smaller ID to conform to USP XXI requirements.

Part Number	Tube ID (mm)	Wall (mm)	Length (mm)	Vial Qty
CG-1841-01	1.5-1.8	0.20	90	100
CG-1841-02	0.8-1.1	0.25	90	100


CG-1182 TLC SPOTTING CAPILLARY TUBES

Non-calibrated tubes made from borosilicate glass for use with aqueous and organic solutions. Available in both 4" and 6" lengths and packed in half pound packages.

Part Number	Length (in)	Approx. Pcs. /Pkg.
CG-1182-01	4	2000
CG-1182-02	6	1300
CG-1182-21	8	1300



CG-1840
CG-1841



CG-1182