


CG-1904 PLATE GLASS, SQUARE, BOROSILICATE

Pyrex® or Duran® borosilicate plate glass. The 1/16", 1/8" and 1/4" thick plates have edges that are cut and swiped. Sizes other than those listed are available. Please call for details.

| Part Number | Size (in) | Thickness (in) | Part Number | Size (in) | Thickness (in) |
|-------------|-----------|----------------|-------------|-----------|----------------|
| CG-1904-36 | 2 X 2 | 1/16 | CG-1904-40 | 8 X 8 | 1/16 |
| CG-1904-01 | 2 X 2 | 1/8 | CG-1904-21 | 8 X 8 | 1/8 |
| CG-1904-02 | 2 X 2 | 1/4 | CG-1904-22 | 8 X 8 | 1/4 |
| CG-1906-37 | 3 X 3 | 1/16 | CG-1904-41 | 10 X 10 | 1/16 |
| CG-1904-06 | 3 X 3 | 1/8 | CG-1904-26 | 10 X 10 | 1/8 |
| CG-1904-07 | 3 X 3 | 1/4 | CG-1904-27 | 10 X 10 | 1/4 |
| CG-1904-38 | 4 X 4 | 1/16 | CG-1904-42 | 12 X 12 | 1/16 |
| CG-1904-11 | 4 X 4 | 1/8 | CG-1904-31 | 12 X 12 | 1/8 |
| CG-1904-12 | 4 X 4 | 1/4 | CG-1904-32 | 12 X 12 | 1/4 |
| CG-1904-13 | 4 X 4 | 1/2 | | | |
| CG-1904-39 | 6 X 6 | 1/16 | | | |
| CG-1904-16 | 6 X 6 | 1/8 | | | |
| CG-1904-17 | 6 X 6 | 1/4 | | | |
| CG-1904-18 | 6 X 6 | 1/2 | | | |




CG-1904

CG-1906 PLATE GLASS, DISC, BOROSILICATE

Pyrex® or Duran® borosilicate glass disc. Edges on all sizes are ground and chamfered. Sizes other than those listed are available. Please call for details.

| Part Number | Dia (in) | Thickness (in) | Part Number | Dia (in) | Thickness (in) |
|-------------|----------|----------------|-------------|----------|----------------|
| CG-1906-41 | 1 | 1/16" | CG-1906-45 | 6 | 1/16" |
| CG-1906-01 | 1 | 1/8" | CG-1906-21 | 6 | 1/8" |
| CG-1906-02 | 1 | 1/4" | CG-1906-22 | 6 | 1/4" |
| CG-1906-42 | 2 | 1/16" | CG-1906-23 | 6 | 1/2" |
| CG-1906-06 | 2 | 1/8" | CG-1906-46 | 8 | 1/16" |
| CG-1906-07 | 2 | 1/4" | CG-1906-26 | 8 | 1/8" |
| CG-1906-08 | 2 | 1/2" | CG-1906-27 | 8 | 1/4" |
| CG-1906-43 | 3 | 1/16" | CG-1906-47 | 10 | 1/16" |
| CG-1906-11 | 3 | 1/8" | CG-1906-31 | 10 | 1/8" |
| CG-1906-12 | 3 | 1/4" | CG-1906-32 | 10 | 1/4" |
| CG-1906-13 | 3 | 1/2" | CG-1906-48 | 12 | 1/16" |
| CG-1906-44 | 4 | 1/16" | CG-1906-36 | 12 | 1/8" |
| CG-1906-16 | 4 | 1/8" | CG-1906-37 | 12 | 1/4" |
| CG-1906-17 | 4 | 1/4" | | | |
| CG-1906-18 | 4 | 1/2" | | | |




CG-1906

CG-1179 PLATE GLASS CUTTER, DIAMOND TIP

Nickel plated brass cutter has a diamond tip with a red wooden handle. Rack has four slots that will accommodate plates from 1/16" to 13/16" (2 to 4.5mm) thick.

| Part Number | Description |
|-------------|------------------------------|
| CG-1179-20 | Glass Cutter, Diamond Tipped |




CG-1179

QUARTZ PRODUCTS

CGQ-0310 GRADED SEAL, QUARTZ TO BOROSILICATE

Graded seals have 3" of quartz sealed to 4" of borosilicate.

| Part Number | ID x OD (mm) | Part Number | ID x OD (mm) |
|-------------|--------------|-------------|--------------|
| CGQ-0310-01 | 2 x 4 | CGQ-0310-11 | 12 x 14 |
| CGQ-0310-02 | 3 x 5 | CGQ-0310-12 | 13 x 15 |
| CGQ-0310-03 | 4 x 6 | CGQ-0310-13 | 14 x 16 |
| CGQ-0310-04 | 5 x 7 | CGQ-0310-14 | 15 x 17 |
| CGQ-0310-05 | 6 x 8 | CGQ-0310-15 | 16 x 18 |
| CGQ-0310-06 | 7 x 9 | CGQ-0310-16 | 17 x 19 |
| CGQ-0310-07 | 8 x 10 | CGQ-0310-17 | 18 x 20 |
| CGQ-0310-08 | 9 x 11 | CGQ-0310-18 | 20 x 22 |
| CGQ-0310-09 | 10 x 12 | CGQ-0310-20 | 22 x 25 |
| CGQ-0310-10 | 11 x 13 | CGQ-0310-21 | 25 x 28 |
| | | CGQ-0310-22 | 27 x 30 |




CGQ-0310

CGQ-0425 STOPCOCK, QUARTZ, STRAIGHT BORE, TEFLON PLUG

Straight bore stopcock made from clear fused quartz. For replacement Teflon plugs see CG-424.

| Part Number | Bore Size (mm) | Plug Size | Sidearm OD (mm) |
|-------------|----------------|-----------|-----------------|
| CGQ-0425-01 | 2 | 11/25 | 8 |
| CGQ-0425-02 | 4 | 15.2/30 | 10 |




CGQ-0425

CGQ-0600 PLATE GLASS, DISC, QUARTZ


Commercial quality ground and polished clear fused quartz discs. Edges have a ground finished on all sizes. **SIZES OTHER THAN THOSE LISTED ARE AVAILABLE. PLEASE CALL FOR DETAILS.**

| Part Number | Dia (in) | Thickness (in) | Part Number | Dia (in) | Thickness (in) |
|-------------|----------|----------------|-------------|----------|----------------|
| CGQ-0600-01 | 1 | 1/16 | CGQ-0600-13 | 4 | 1/16 |
| CGQ-0600-02 | 1 | 1/8 | CGQ-0600-14 | 4 | 1/8 |
| CGQ-0600-03 | 1 | 1/4 | CGQ-0600-15 | 4 | 1/4 |
| CGQ-0600-04 | 1.5 | 1/16 | CGQ-0600-16 | 5 | 1/16 |
| CGQ-0600-05 | 1.5 | 1/8 | CGQ-0600-17 | 5 | 1/8 |
| CGQ-0600-06 | 1.5 | 1/4 | CGQ-0600-18 | 5 | 1/4 |
| CGQ-0600-07 | 2 | 1/16 | CGQ-0600-19 | 6 | 1/16 |
| CGQ-0600-08 | 2 | 1/8 | CGQ-0600-20 | 6 | 1/8 |
| CGQ-0600-09 | 2 | 1/4 | CGQ-0600-21 | 6 | 1/4 |
| CGQ-0600-10 | 3 | 1/16 | | | |
| CGQ-0600-11 | 3 | 1/8 | | | |
| CGQ-0600-12 | 3 | 1/4 | | | |




CGQ-0600


QUARTZ PRODUCTS



CGQ-0620



CGQ-0640



CGQ-0660



CGQ-0690

CGQ-0620 PLATE GLASS, SQUARE, QUARTZ

Commercial quality ground and polished clear fused quartz plate. Edges have a ground finished on all sizes. **SIZES OTHER THAN THOSE LISTED ARE AVAILABLE. PLEASE CALL FOR DETAILS.**

| Part Number | Size (in) | Thickness (in) | Part Number | Size (in) | Thickness (in) |
|-------------|-----------|----------------|-------------|-----------|----------------|
| CGQ-0620-01 | 1 X 1 | 1/16 | CGQ-0620-13 | 4 X 4 | 1/16 |
| CGQ-0620-02 | 1 X 1 | 1/8 | CGQ-0620-14 | 4 X 4 | 1/8 |
| CGQ-0620-03 | 1 X 1 | 1/4 | CGQ-0620-15 | 4 X 4 | 1/4 |
| CGQ-0620-04 | 1.5 X 1.5 | 1/16 | CGQ-0620-16 | 5 X 5 | 1/16 |
| CGQ-0620-05 | 1.5 X 1.5 | 1/8 | CGQ-0620-17 | 5 X 5 | 1/8 |
| CGQ-0620-06 | 1.5 X 1.5 | 1/4 | CGQ-0620-18 | 5 X 5 | 1/4 |
| CGQ-0620-07 | 2 X 2 | 1/16 | CGQ-0620-19 | 6 X 6 | 1/16 |
| CGQ-0620-08 | 2 X 2 | 1/8 | CGQ-0620-20 | 6 X 6 | 1/8 |
| CGQ-0620-09 | 2 X 2 | 1/4 | CGQ-0620-21 | 6 X 6 | 1/4 |
| CGQ-0620-10 | 3 X 3 | 1/16 | | | |
| CGQ-0620-11 | 3 X 3 | 1/8 | | | |
| CGQ-0620-12 | 3 X 3 | 1/4 | | | |

CGQ-0640 MICROSCOPE SLIDES, QUARTZ

Ground and polished clear quartz microscope slides. All sizes are 1mm thickness.

| Part Number | Dimensions (mm) / (in.) |
|-------------|-------------------------|
| CGQ-0640-01 | 75 x 25 (3" x 1") |
| CGQ-0640-02 | 50 x 25 (2" x 1") |
| CGQ-0640-03 | 25 x 25 (1" x 1") |

CGQ-0660 COVER SLIPS, QUARTZ

Ground and polished clear quartz microscope slide cover slips are ultra-thin (0.25mm) thickness.

| Part Number | Dimensions (mm) |
|-------------|-----------------|
| CGQ-0660-01 | 19 x 19 |

CGQ-0690 QUARTZ WOOL

Useful as high temperature packing in distillation columns or as a bed support for high purity liquid chromatography. Available in fine or coarse porosity. **SOLD IN 1/2lb PACKAGES ONLY.**

| Part Number | Description | Fiber Size |
|-------------|-------------|--------------|
| CGQ-0690-01 | Fine | 4-8 Microns |
| CGQ-0690-02 | Coarse | 8-15 Microns |