

SHAKER



CG-1980



CG-1980-25



CG-1980-26



CG-1980-27

CG-1980 IKA® MINI SHAKER, DIGITAL

Compact, universal small shaker suitable for shaking tasks with all small vessels and microtiter plates.

Supplied complete with the following attachments: (1ea) CG-1980-25 Universal Attachment, (1ea) CG-1980-26 Universal Attachment for Foam Inserts, (1ea) CG-1980-27 Microtiter Plate Attachment, (1ea) CG-1980-29 One-Hand Insert, and (1ea) CG-1980-30 6 Test Tube Insert.

Features:

- Wide Range Of Attachments
- Attachment Detection
- Continuous Or Touch Operation (With Standard Attachment)
- Two Operating Modes: **Mode A:** (Safe Mode With Attachment Detection). The Maximum Speed Of 3000 RPM Is Only Reached With The Standard Attachment In Touch Mode. When Using Other Attachments, The Speed Is Limited To 1300 RPM. **Mode B:** (Without Attachment Detection). A Speed Range Of 3000 RPM Is Possible With All Attachments.
- Stable In All Speed Ranges
- Sturdy Zinc Die Cast Casing

Specifications:

Shaking Movement:	Orbital
MAX Permitted Shaking Weight with Attachment (lbs):	1.1
Speed Range (RPM):	0 to 3000
Speed Display	Digital
Dimensions W x D x H (mm):	148 x 205 x 63
Weight (lbs):	6.38

Part Number	Description
CG-1980-11	Mini Shaker, Complete

Replacement Items Included with Shaker:

Part Number	Description
CG-1980-25	Standard Attachment for Test Tubes And Small Vessels
CG-1980-26	Universal Attachment for Foam Inserts
CG-1980-27	Microtiter Plate Attachment
CG-1980-29	One-Hand Insert, Foam
CG-1980-30	6 Test Tube Insert, 12mm Diameter, Foam

Optional Items for Shaker:

Part Number	Description
CG-1980-33	14 Test Tube Insert, 10mm Diameter, Foam
CG-1980-34	4 Test Tube Insert, 16mm Diameter, Foam
CG-1980-35	Blank Insert for Adding Any Number of Holes, Foam