

SLEEVES

CG-116 SLEEVE, TEFLON


Reusable 0.005" thick Teflon sleeve for standard taper joints provides a grease-free seal.

MEDIUM LENGTH

Part Number	Joint Size
CG-116-02	10/18
CG-116-04	14/20
CG-116-05	19/22
CG-116-06	24/25
CG-116-07	29/26
CG-116-09	40/35

FULL LENGTH

Part Number	Joint Size
CG-116-11	10/30
CG-116-13	14/35
CG-116-14	19/38
CG-116-15	24/40
CG-116-16	29/42
CG-116-17	34/45
CG-116-18	40/50
CG-116-19	45/50
CG-116-21	55/50
CG-116-22	71/60



CG-116

CG-136 SLEEVE, TEFLON, RIBBED, GRIPPING RING

Reusable 0.008" thick Teflon sleeve for standard taper joints has a knurled, reinforced gripping ring for easy removal. Sleeve has reinforcement ribs that provide a vacuum tight seal to 1 mbar and also prevent the sleeve from sticking.

Part Number	Joint Size
CG-136-01	14/20
CG-136-02	19/22
CG-136-03	24/40
CG-136-04	29/42
CG-136-05	34/45
CG-136-06	45/50
CG-136-07	55/50




CG-136

CG-137 SLEEVE, TEFLON, FLURO-KEM

A 0.002" thick Teflon sleeve for standard taper joints provides a grease-free seal. **SOLD 3/PKG.**

Part Number	Joint Size
CG-137-02	14/35
CG-137-03	19/38
CG-137-04	24/40
CG-137-05	29/42
CG-137-06	34/45
CG-137-07	45/50



CG-137

CG-110 SPRING, STAINLESS STEEL

Stainless steel spring used to connect apparatus having hooks. All springs have an enlarged 1/4" ID loop at each end for easy connection to almost any size hook. Diameter of the coils on all sizes is 0.190". Length listed is the length between the loops. **SOLD 12/PKG.**

Part Number	Length of Coil (mm)
CG-110-01	6.35
CG-110-02	18.0
CG-110-03	25.0
CG-110-04	32.0



CG-110

Visit Our Web Site: www.gpelimited.co.uk