


THERMOMETERS / TUBES



CG-3508-N

CG-3508-N THERMOMETER, NON-MERCURY, 10/18 JOINT

Thermometer is mineral spirit filled to eliminate workplace and environmental hazards associated with mercury. All ranges have 1° subdivisions. Thermometer has a white back and a 6 to 7mm outer diameter with a standard taper 10/18 inner ground joint. Immersion is measured from the bottom of the joint to the bottom of the mercury bulb.


Part Number	Immersion (mm)	Range (°C)
CG-3508-N-01	25	-20 to 150
CG-3508-N-02	25	-10 to 250
CG-3508-N-03	51	-20 to 150
CG-3508-N-04	51	-10 to 250
CG-3508-N-06	76	-10 to 250

CG-4013

CG-4013 TUBE, CENTRIFUGE, 50mL, UNGRADUATED

Centrifuge tube with a short, conical bottom and 24-415 GPI thread. Fits I.E.C. shield number 320. Supplied complete with a Poly-Seal® polypropylene lined phenolic screw cap.

Part Number	Capacity (mL)	GPI Thread Size	Approx OD x Length (mm)
CG-4013-01	50	24-415	28 x 155




CG-4022

CG-4022 TUBE, CONNECTING, T-SHAPE

T-shaped tube having serrated ends for connection of rubber, latex or Tygon® tubing.

Part Number	Tubing OD (mm)	Hose Connection
CG-4022-01	8.0	CG-300-01
CG-4022-02	10.0	CG-300-02
CG-4022-03	12.7	CG-300-03



CG-4023

CG-4023 TUBE, CONNECTING, Y-SHAPE

Y-shaped tube having serrated ends for connection of rubber, latex or Tygon® tubing.

Part Number	Tubing OD (mm)	Hose Connection
CG-4023-01	8.0	CG-300-01
CG-4023-02	10.0	CG-300-02
CG-4023-03	12.7	CG-300-03