

External pressure pump. Bearing out of the liquid circuit. All parts in the liquid circuit and housing made of stainless steel. Operation either manual (on/off) or external via 3-line-connector. Cable (Option) for the potential free contact of Unistat Tango and the big power Unistats, all Compatible Control Thermostats and Unichillers.

Technical data according to DIN 12876

Operating temperature range	-120...300 °C
Circulation pump	Yes
max. delivery of pressure pump	70 l/min
max. delivery pressure	2.5 bar
Pump connection	M30x1,5 male
Overall dimensions WxDxH **	215x457x295 mm
Net weight	26 kg
sound pressure level	64 dB(A)
Power supply (3 Phase)	400V 3~N 50Hz
max. current (3 Phase)	2,4 A
Fuse (3 phase)	3x4 A
min. ambient temperature	5 °C
max. ambient temperature	32 °C

**Order-No.: 527.0020****from Serial-No.:****146630****1.2/13**

Technical details and dimensions are subject to change. No liability is accepted for errors or omissions.

Output data valid for: Room temperature 20°C

** Please respect space requirements.