

External pressure pump. Bearing out of the liquid circuit. All parts in the liquid circuit and housing made of stainless steel. Operation either manual (on/off) or external via 3-line-connector. Cable (Option) for the potential free contact of Unistat Tango and the big power Unistats, all Compatible Control Thermostats and Unichillers.

Technical data according to DIN 12876

Operating temperature range	-120...300 °C	
Circulation pump	Yes	
max. delivery of pressure pump	100 l/min	Order-No.: 527.0016
max. delivery pressure	3,7 bar	
Delivery at 0,5 bar	99 l/min	
Delivery at 1,0 bar	98 l/min	
Delivery at 1,5 bar	97 l/min	
Delivery at 2,0 bar	96 l/min	
Delivery at 2,5 bar	92 l/min	
Delivery at 3,0 bar	76 l/min	
Delivery at 3,5 bar	45 l/min	
Pump connection	M30x1,5 male	
Overall dimensions WxDxH **	215x457x295 mm	
Net weight	26 kg	
Power supply (3 Phase)	460V 3~ 60Hz	
max. current (3 Phase)	3 A	
Fuse (3 phase)	3x6 A	
min. ambient temperature	5 °C	
max. ambient temperature	32 °C	

from Serial-No.: **170704** **1.2/15**

Technical details and dimensions are subject to change. No liability is accepted for errors or omissions.

Output data valid for: Room temperature 20°C

** Please respect space requirements.