

External pressure pump. Bearing out of the liquid circuit. All parts in the liquid circuit and housing made of stainless steel. Operation either manual (on/off) or external via 3-line-connector. Cable (Option) for the potential free contact of Unistat Tango and the big power Unistats, all Compatible Control Thermostats and Unichillers.

**Technical data according to DIN 12876**

---

Operating temperature range	-120...300 °C	
Circulation pump	Yes	
max. delivery of pressure pump	70 l/min	<b>Order-No.: 527.0015</b>
max. delivery pressure	1.5 bar	
Pump connection	M30x1,5 male	
Overall dimensions WxDxH **	215x457x295 mm	
Net weight	20 kg	
Power supply (3 Phase)	460V 3~ 60Hz	
max. current (3 Phase)	2 A	
Fuse (3 phase)	3x4 A	

---

**from Serial-No.:****180263****1.0/13**

---

Technical details and dimensions are subject to change. No liability is accepted for errors or omissions.