

External pressure pump. Bearing out of the liquid circuit. All parts in the liquid circuit and housing made of stainless steel. Operation either manual (on/off) or external via 3-line-connector. Cable (Option) for the potential free contact of Unistat Tango and the big power Unistats, all Compatible Control Thermostats and Unichillers.

Technical data according to DIN 12876

Operating temperature range	-120...300 °C	
Circulation pump	Yes	
max. delivery of pressure pump	70 l/min	Order-No.: 527.0018
max. delivery pressure	2.5 bar	
Pump connection	M38x1,5 AG	
Overall dimensions WxDxH **	457x201x320 mm	
Power supply (3 Phase)	460V 3~ 60Hz	
min. ambient temperature	5 °C	
max. ambient temperature	32 °C	

from Serial-No.:**1.0/11**

Technical details and dimensions are subject to change. No liability is accepted for errors or omissions.

Output data valid for: Room temperature 20°C

** Please respect space requirements.