

External pressure pump. Bearing out of the liquid circuit. All parts in the liquid circuit and housing made of stainless steel. Operation either manual (on/off) or external via 3-line-connector. Cable (Option) for the potential free contact of Unistat Tango and the big power Unistats, all Compatible Control Thermostats and Unichillers.

### Technical data according to DIN 12876

Operating temperature range	-120...300 °C
max. delivery of pressure pump	70 l/min
max. delivery pressure	1.5 bar
Pump connection	M38x1,5 AG
Overall dimensions WxDxH **	215x457x295 mm
Net weight	25 kg
Power supply (3 Phase)	460V 3~ 60Hz
max. current (3 Phase)	2,8 A
Fuse (3 phase)	3x4 A
min. ambient temperature	5 °C
max. ambient temperature	32 °C



**Order-No.: 527.0017**

**from Serial-No.:**

**178770**

**1.0/13**

Technical details and dimensions are subject to change. No liability is accepted for errors or omissions.

Output data valid for: Room temperature 20°C

\*\* Please respect space requirements.