

Save Time



Cost Effective



Easy-to-use



AUTO Sample - 60



GPE Scientific Limited

**Greaves Way Industrial Estate
Stanbridge Road, Leighton Buzzard
Bedfordshire, LU7 4UB
United Kingdom**

**Tel: 44 (0) 1525 382277
Fax: 44 (0) 1525 382263**

**Email: sales@gpescientific.co.uk
Website: www.gpescientific.co.uk**



nanalysis

TECHNICAL SPECIFICATIONS

**System Requirements: NMReady-60e or 60PRO
Dimensions: 10.3 x 10.3 x 14.2" / 26 x 26 x 36 cm
Weight: 12 lbs/5.5 kg
Voltage: 24 V
Time per Tube Change: ~10 sec
Carousels: 2**



Maximize your efficiency by automating routine work.

Introducing the affordable, easy-to-use AUTOsample-60.



The AUTOsample-60 mounts on top of any NMReady-60e or 60PRO for painless interconversion between autosampler and manual mode, while still maintaining the notably minimal footprint. Data acquisition and processing is controlled at the touch of a finger through the onboard software.

EASY CONFIGURATION

Simple mounting on top of NMReady-60

Interchangeable 25 NMR tube carousels (standard, constricted or J-Young tubes)

Operated through onboard software via NMReady USB connectors

Robust and reliable tube changes ensure safety and accuracy

Optional sample tempering

SIMPLE OPERATION

Use onboard software to configure data acquisition of up to 25 different samples at one time.

Set up structure elucidation or conformation protocols per tube

e.g., 1H, COSY, HSQC, 13 C{1H}

QA/QC for series of samples

